

QUICK LOOK

OAE + **ABR**

PRODUCT Guide



Otoport
ADVANCE

Otoport
CLINICAL

Otocheck
DP

Echoport 292

Otocheck
TB

Otoport
TB SCREENER

Otocheck
ABR

Anniversary
Trade in special offer

Quick Look Product Guide



Echoport 292

(Laptop or Desktop PC required)



Otoport ADVANCE



Otoport DP CLINICAL



Otocheck DP



Otocheck TE



Otoport TE SCREENER



Otocheck ABR

	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
--	---------	----------	----------	----------	----------	----------	----------

USE	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
NEWBORN SCREENING	✓	✓	✓	✓	✓	✓	✓
PEDIATRIC ASSESSMENT	✓	✓	✓	✓	✓	✓	✓
CLINICAL AUDIOLOGY	✓	✓	✓	✓	✓	✓	✓
RESEARCH	✓	✓	✓	✓	✓	✓	✓

DP OAE	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
TEST CONFIGURABILITY	4 Protocols, 14 Parameters	7 Protocols, 16 Parameters	6 Protocols, 11 Parameters	1 Protocol, 7 Parameters	-	-	-
FREQUENCY RANGE	0.8 to 10.5 kHz	1 to 8 kHz, 13 Frequencies	1 to 8 kHz, 12 Frequencies	2, 3, 4, 6 kHz	-	-	-
FREQ. RESOLUTION	1,2,3,4,8 pts/octave	Up to 16 pts/octave	Up to 4 pts/octave	1/2 octave	-	-	-
FREQUENCY RATIO	1.04 to 1.5 (Δ0.02)	-	-	-	-	-	-
STIMULUS LEVELS (dB SPL)	40 to 80 (Δ1dB)	40 to 75 (Δ5dB)	60/50 to 70/70 (6 Presets)	60/50 to 70/70 (6 Presets)	-	-	-
PASS CRITERIA	10 Parameters	9 Parameters	3 Parameters	2, 3, 4 Bands	-	-	-
PASS SNR	1 to 14dB per 1/2 octave (Δ1dB)	1 to 14dB per band (Δ1dB)	1 to 10dB per band (Δ1dB)	1 to 10dB per band (Δ1dB)	-	-	-

TE OAE	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
TEST CONFIGURABILITY	4 Protocols, 10 Parameters	3 Protocols, 16 Parameters	2 Protocols, 10 Parameters	-	1 Protocol, 7 Parameters	3 Protocols, 16 Parameters	-
FREQUENCY RANGE	1 to 6 kHz, 100 pts	1, 1.5, 3, 4, 5, 6 kHz	1, 1.5, 2, 3, 4kHz	-	1, 1.5, 2, 3, 4kHz	1, 1.5, 3, 4, 5, 6kHz	-
FREQ. RESOLUTION	50Hz, 80Hz or 1/2 octave	80Hz or 1/2 Octave bandwidth	1/2 Octave bandwidth	-	1/2 Octave bandwidth	1/2 Octave bandwidth	-
STIMULUS LEVELS (dB SPL _{pe})	40 to 90 (Δ1dB)	60 to 90 (Δ1db)	60 to 84 (Δ1db)	-	60 to 84 (Δ1dB)	60 to 90 (Δ1dB)	-
PASS CRITERIA	10 Parameters	7 Parameters	2, 3, or 4 1/2 octave, wide or narrow bands to pass	-	2, 3, or 4 1/2 octave, wide or narrow bands to pass	7 Parameters	-
PASS SNR	1 to 14dB per 1/2 octave (Δ1dB)	1 to 14dB (Δ1dB)	6dB (Fixed)	-	6dB (Fixed)	1 to 14dB (Δ1dB)	-

ABR	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
TEST CONFIGURABILITY	-	OPTIONAL ABR	OPTIONAL ABR	OPTIONAL ABR	OPTIONAL ABR	OPTIONAL ABR	2 Protocols, 9 Parameters
STIMULUS LEVELS (dBnHL)	-	-	-	-	-	-	5 to 60dBnHL (Δ1dB)
STIMULUS TYPE	-	-	-	-	-	-	Chirp and Click

HARDWARE	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
BATTERY LIFE	No Batteries required	250 tests per battery charge	250 tests per battery charge	250 tests per battery charge	250 tests per battery charge	250 tests per battery charge	250 tests per battery charge
DATABASE STORAGE	Unlimited	5000 tests, 1000 patients	1000 tests, 1000 patients	500 tests, 500 patients	500 tests, 500 patients	5000 tests, 1000 patients	500 tests, 500 patients
LANGUAGES	English, French, German, Italian, Portuguese	English	English, Spanish, French, Portuguese, Italian, German, Turkish, Polish, Hungarian	English, Spanish, French, Portuguese, Italian, German, Turkish, Polish, Hungarian	English, Spanish, French, Portuguese, Italian, German, Turkish, Polish, Hungarian	English	English, Spanish, French, Portuguese, Italian, German, Turkish, Polish, Hungarian
PATIENT DEMOGRAPHIC <i>Further patient demographic fields available via Otolink</i>	Extensive patient database fields	ID, name, DOB, gender, mother, notes, risk, location, facility, patient status, consent	ID/notes, name, DOB, gender, location, facility, patient status	ID/notes, name, DOB, gender, location, facility, patient status	ID/notes, name, DOB, gender, location, facility, patient status	ID, name, DOB, gender, mother, notes, risk, location, facility, patient status, consent	ID/notes, name, DOB, gender, location, facility, patient status
CPT codes	92588, 92587, 92558	92588, 92587, 92558 (92586 with ABR option)	92588, 92587, 92558 (92586 with ABR option)	92587, 92558 (92586 with ABR option)	92558 (92586 with ABR option)	92558 (92586 with ABR option)	92586 (92558, 92587, 92588 with OAE option)

KEY FEATURES	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
<p>* 2 YEARS INSTRUMENT WARRANTY</p> <p>** 3 YEARS INSTRUMENT WARRANTY</p> <ul style="list-style-type: none"> • SOAE synchronous spontaneous emission search • DPgram high resolution 8pts/oct (27 frequencies) • DP Growth I/O with scissor mode • DP&TE binaural recording*** • TE Contralateral suppression*** • TE Quickscreen, non-linear, linear, tone pip 	<ul style="list-style-type: none"> • DP Zoom to 16pts/oct • DP Growth I/O with scissor mode • Smart SNR, configurable -15 to 0 • Configurable individual pass criteria • Multi test macro (4xprotocols) • Extensive patient database fields • Dedicated QA test area & history view • TEOAE Waveform • Plus Otoport DP Clinical features 	<ul style="list-style-type: none"> • DPgram line & bar graph display • Data Table display • One key select Clinical/Screening protocol (unique to DP Clinical) • Linked test macro function • TEOAE Waveform (Add Option) • Pediatric distraction display • Plus Otocheck features 	<ul style="list-style-type: none"> • Rapid DPOAE screening mode • Quick save mode 2 key press test • Individual frequency progress display • Easy-to-use pass/refer result • Create multiple users with Admin management and security • Programmable worklist • Probe test function and ID tracking 	<ul style="list-style-type: none"> • Gold standard TEOAE Quickscreen • Quick save mode 2 key press test • Choice in progress display • Easy-to-use pass/refer result • Create multiple users with Admin management and security • Programmable worklist • Probe test function and ID tracking 	<ul style="list-style-type: none"> • Gold standard TEOAE Quickscreen • Configurable individual pass criteria • Fast pass/refer result • Extensive patient database fields • Create multiple users with Admin management and security • Programmable worklist • Dedicated QA test area & history view 	<ul style="list-style-type: none"> • Superior pretest environment, with powerline and impedance analysis • Clear test progression notification • Fast pass/refer result • Choice in stim delivery ear cup/inset • Real time ABR waveform view • Create multiple users with Admin management and security • Programmable worklist 	

ADDITIONAL OPTIONS	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
Add Pediatric Distraction Display	-	INCLUDED	INCLUDED	INCLUDED AS SPECIAL OFFER	INCLUDED AS SPECIAL OFFER	-	-
Add Wireless Module	-	INCLUDED	INCLUDED AS SPECIAL OFFER	INCLUDED AS SPECIAL OFFER	INCLUDED AS SPECIAL OFFER	INCLUDED AS SPECIAL OFFER	INCLUDED AS STANDARD
Add Docking Station	-	INCLUDED	INCLUDED AS SPECIAL OFFER	OPTION	OPTION	OPTION	-
Add Mini Printer	-	OPTION	OPTION	OPTION	OPTION	OPTION	OPTION
Add More	-	ADD Noah	ADD TEOAE, ADD ILOV6, ADD Noah	ADD Noah	ADD Noah	ADD Noah	ADD Binaural testing, ADD Noah

TRADE IN SPECIAL	DESKTOP	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD	HANDHELD
STANDARD LIST PRICE	FROM \$7,499	FROM \$7,295	FROM \$4,999	FROM \$4,195	FROM \$4,195	FROM \$4,490	FROM \$9,500
Otodynamics legacy product	\$1,000 off list price	\$1,000 off list price	\$1,000 off list price	\$750 off list price	\$750 off list price	\$750 off list price	\$1,000 off list price
Other manufacturer legacy products	\$750 off list price	\$750 off list price	\$750 off list price	\$500 off list price	\$500 off list price	\$500 off list price	\$750 off list price



ANNIVERSARY YEAR

Longevity equals experience

Founded in 1988 by the discoverer of Otoacoustic emissions, Dr David T Kemp FRS. Otodynamics was the first OAE company. Being a privately owned and managed company it's free to give top priority to research and development. Its in-house engineers and scientists are experts dedicated to quality, reliability and performance. Otodynamics stands out as the gold standard, consistently the first choice in Audiology.



To celebrate 30 years of Otodynamics
See our 4 special offers available 1st April - 31st August

1 Trade in*
UP TO **\$1,000** Off List Price
*T&C apply

2 Free pediatric distraction display with all Otocheck orders
worth **\$250**

3 Free Docking station upgrade with Otoport DP clinical orders
worth **\$270**

4 Free Wireless Module with All Otoport/Otocheck orders
worth **\$490**



For further information, please contact us:
www.otodynamics.com
info@otodynamics.com

* Otoport Instruments ** Echoport Instruments *** Echoport ILO292 II



ANNIVERSARY EDITION



Otodynamics Ltd.

30-38 Beaconsfield Road Hatfield Herts AL10 8BB UK

t. 1-800-659-7776 e. info@otodynamics.com

www.otodynamics.com